


Lab name: Canvas Labs  
Client: [REDACTED]  
Analysis date: 11/23/2017 16:41:47  
Method: Syringe Injection  
Column: RESTEK 15METER MXT-1  
Sample: Sovrin Cart SLH (Distillate)

Temperature program:

Int temp	Hold	Ramp	Final temp
140.00	0.000	20.000	260.00
100.00	0.000	20.000	140.00

Events:

Time	Event
0.000	ZERO
0.000	SOUND
2.000	INTEG IMMEDIATE



Component Internal Units

Internal Standard	100.0
CBD	0.0%
d-8	5.2%
THC	85.1%
CBN	5.4%
Total	95.6%